



Eppendorf – In touch with life[®]



Eppendorf Safe-Lock Tubes® BioBased 0.5/1.5/2.0 mL

Hamburg, June 2025

Public Use

Eppendorf SE

Agenda


1 General Information	<u>04-07</u>
2 Customers to Address	<u>08</u>
3 What Covers an ESG Report	<u>09</u>
4 Customer Voice	<u>10-11</u>
5 Marketing Material	<u>12</u>
6 Order Numbers	<u>13</u>

General Information

- Gate 5 Date: June 13th 2025

A large, light beige arrow pointing to the right, containing the text "Gate 5: June 13th".

Gate 5: June 13th

A large, light beige arrow pointing to the right, containing the text "Go Live EMEA: June 23rd".

Go Live EMEA:
June 23rd

General Information

- 100 % BioBased, made from 100 % renewable resources (2nd generation feedstock such as waste cooking oil)
- 0.5 mL, 1.5 mL and 2.0 mL format
- colorless
- Available in Eppendorf Quality, and PCR clean
- Performance is same as the fossil-based Safe-lock Tubes (according to individual feedback of customers during customer test)



General Information

- Centrifugation up to 30.000 x g
- Temperature range – 86 °C to 100 °C
- Autoclavable with open lid (120 °C, 20 min)
- Safe incubation and storage (no opening) due to Safe-Lock lid
- No adaption of assays necessary when switching from Eppendorf Safe-Lock Tubes® to Biobased Tubes, as they perform similar



General Information

Packaging

- Eppendorf Safe-Lock Tubes® BioBased are packed in resealable bags, 500 tubes, colorless
- For better identification, an explanatory additional ISCC label will be placed next to the product label on the resealable bags
- The accompanying text is formulated according to the ISCC guideline and must not be changed



Customers to Address

Companies publishing an ESG report

An ESG report shows the environmental, social and governance measures of companies to create a better future by reducing e.g. use of fossil resources. Companies publishing this report have a general interest and pressure to reduce their carbon footprint in all levels of their business.

What covers an ESG report?

How a companies' activities impact the environment, society and employees on different aspects such as

- > Energy consumption
- > CO₂ emissions
- > Compliance with human rights
- > Working conditions
- > Corporate ethics

Mandatory for (European) Companies with a

- > Balance sheet total over 25 million euros
- > Net sales higher than 50 million euros
- > More than 250 employees

- > Also, Eppendorf, Fisher Scientific, etc. must publish an ESG report

Customer Voices

"Switching to bio-based consumables is a big step toward making lab plastics more sustainable – without sacrificing quality or performance. With Eppendorf BioBased Tubes, we've reduced our environmental impact while maintaining the reliability we trust. It's an easy and effective move toward a greener future in the lab."

- José Pedro Soares Maio (University of Namur/URBM: Unité de Recherche en Biologie des Micro-organismes)

Customer Voices

"As a scientist passionate about innovation and sustainability, I was excited to put Eppendorf BioBased Tubes to the test in our lab. After thorough evaluation, we found that they perform just as well as conventional fossil-based Eppendorf Tubes across all our applications: molecular biology, cloning, centrifugation, heat applications."

- José Pedro Soares Maio (University of Namur/URBM: Unité de Recherche en Biologie des Micro-organismes)

Marketing Material

- > Images Landing Page www.Eppendorf.com/biobased
- > Social Media Image + 3 Text Suggestions
- > Newsletter Text
- > Customer Presentation
- > AppNote 495 and 496 “Costs of reproducibility” + corresponding Sales Guide 1-Pager

Order Numbers

Description	Order No.	Reference (Fossil-based Eppendorf Safe-Lock Tubes)	
Safe-Lock Tubes® 0.5 mL, BioBased, Eppendorf Quality, 500 tubes	0030123042	0030121023	
Safe-Lock Tubes® 1.5 mL, BioBased, Eppendorf Quality, 1.000 tubes	0030123085	0030120086	
Safe-Lock Tubes® 2.0 mL, BioBased, Eppendorf Quality, 1.000 tubes	0030123166	0030120094	
Safe-Lock Tubes® 0.5 mL, BioBased, PCR Clean, 500 tubes	0030123069	0030123301	
Safe-Lock Tubes® 1.5 mL, BioBased, PCR Clean, 1.000 tubes	0030123140	0030123328	
Safe-Lock Tubes® 2.0 mL, BioBased, PCR Clean, 1.000 tubes	0030123182	0030123344	
Safe-Lock Tubes® 0.5 mL, BioBased, Eppendorf Quality, 500 tubes	0030123050	022363611	US Version
Safe-Lock Tubes® 1.5 mL, BioBased, Eppendorf Quality, 500 tubes	0030123093	022363204	US Version
Safe-Lock Tubes® 2.0 mL, BioBased, Eppendorf Quality, 500 tubes	0030123174	022363352	US Version
Safe-Lock Tubes® 0.5 mL, BioBased, PCR Clean, 500 tubes	0030123077	022363719	US Version
Safe-Lock Tubes® 1.5 mL, BioBased, PCR Clean, 500 tubes	0030123158	022363212	US Version
Safe-Lock Tubes® 2.0 mL, BioBased, PCR Clean, 500 tubes	0030123190	022363344	US Version

Go For It





Eppendorf – In touch with life[®]

Thank you

Disclaimer

Eppendorf®, the Eppendorf logo, In touch with life are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright © 2025 by Eppendorf SE.