



Flex up Your Team

Tired of time-consuming, endless pipetting steps?

The epMotion® 96 Flex's design enables fast and easy liquid handling in both 96-well and 384-well plates. Boost your productivity, complete tasks quicker and save time for other experiments:

- > Covering 0.5-1000 µL with only two dispensing heads
- > Semi-automated tip loading/ejecting
- > Optical auto-detect function of tips
- > Significant time savings compared to 8-channel pipettes



www.eppendorf.link/flex-up-your-team

Dealer partner logo

epMotion® 96 Flex: Smart features

Status light bar
Visual feedback of device status.

Exchangeable dispensing heads
Maximum flexibility with ISO 23783-2 calibrated dispensing heads for a wide volume range (0.5-1000 µL).

Semi-automated tip loading/ejecting
For quick and reliable exchange of epT.I.P.S® Motion reload tip trays.

Smart Indexing
384-well plates can be pipetted by smart indexing functionality.

Manually lift worktable up and lock in place:
Easy and quick adjustment to individual labware heights and filling levels by using the 2-position slider brake.

Compact design
Small footprint, fits biosafety cabinets and can be moved between labs.

USB port
For data exchange behind front cover.

Optical auto-detect function of tips
Convenient verification of different tip types for maximum security of correct tip type being loaded.

Lifting table
2-position slider for source and destination labware. Supports a variety of SLAS/ANSI – size labware such as micro plates, reservoirs and thermo-adapters.

Height adjustment
Set height and auto trigger functions for semi-automated aspiration and dispensing.

TipHolder for epMotion® 96 Flex

- > For using epT.I.P.S® Motion Reloads
- > One TipHolder delivered with each epMotion® 96 Flex base device
- > Fits only for epMotion® 96 Flex
- > Reusable and autoclavable





Ordering information	
Description	Ordering no.
epMotion® 96 Flex Base Device with 2-positions slider. Base device without dispensing head. With TipHolder for epMotion® 96 Flex for use of epT.I.P.S® Motion Reloads, 100 - 240 V, 50 - 60 Hz	5069 000 519
Dispensing Head 96-channel 0.5-300 µL in hard shell case for epMotion® 96 Flex	5069 000 527
Dispensing Head 96-channel 5-1000 µL in hard shell case for epMotion® 96 Flex	5069 000 535
TipHolder for epMotion® 96 Flex for use with epT.I.P.S® Motion Reloads	5069 075 004

Technical specifications subject to change. Errors and omissions excepted.

If you have any questions or would like to receive a quote, please contact your local distributor: